

# HALUX N30



medicinelaser<sup>®</sup>  
LATAM



## Small but mighty! HALUX N30

The HALUX N30 is a solid examination fixture with uncompromisingly good light and a very attractive price. The illumination field has a diameter of 18 cm and an illuminance of up to 30,000 lux at 0.5 m which makes it perfectly equipped for the versatile daily routine in the clinic. The good colour rendering properties and the stepless dimming function offer high functionality.



HALUX N30	technical data
luminous power *	30,000 Lux @ 0.5 m
colour temperature *	4400 K
lighting field @ 1.0 m *	Ø 18 cm
colour rendering index	Ra > 93
usage	touch panel on luminaire head
class of protection	II
protection rating	IP 20

special features	
HALUX N30	<ul style="list-style-type: none"> <li>• continuously dimmable</li> <li>• maintenance-free</li> <li>• spring arm</li> </ul>
variations	<ul style="list-style-type: none"> <li>• gynecology</li> <li>• ambulance</li> </ul>
fastening	<ul style="list-style-type: none"> <li>• pin universal</li> <li>• roller stand</li> </ul>
accessories	<ul style="list-style-type: none"> <li>• wall mount</li> <li>• roller stand</li> <li>• fastening hardware</li> <li>• universal mount</li> <li>• table clamp</li> </ul>

# HALUX N30 – variations



### Variety of versions

The HALUX N30 is available in different versions. Depending on your requirements, it can be mounted directly on the wall or fixed to a rail mount with a rail clamp. The ambulance version can be mounted on C-rails. The gynecology version is suitable for connection to the gynecology chair.

### Conceivably simple

The 18 cm large illumination field can be easily controlled via the touch panel on the luminaire head. The illuminance can be dimmed steplessly to adapt to the lighting requirements. Cleaning is also simple. This is ensured by the seamless design.

## HALUX N30 P SGH and SGV

pin version | flex arm with foot joint



## HALUX N30 P SH and SV

pin version | flex arm



## HALUX N30 – variations



## HALUX N30 P F1

pin version | spring arm with foot joint



## HALUX N30 P Amb

pin version | head joint | ambulance



## HALUX N30 F

integrated roller stand



# HALUX N50 and N30 – gynecology



## Lighting solutions for gynecology & urology

We support gynecological and urological clinics as well as proctologists with optimal light. Very good colour rendering properties with Ra values > 93 ensure a natural representation of visible pathologies, tissue and vessels. Unobtrusive spotlights on the examination chair provide ideal light for patient examination.

### HALUX N50-3 P SX Gyn

pin version | flex arm | gynecology

- 50,000 Lux @ 0.5 m
- optional: 5-stage dimmable
- adjustable colour temperature (3300/ 3800/ 4400 K)



#### Securely fastened

rail clamp with brake function

### HALUX N30-1 P SGV Gyn

pin version | flex arm | gynecology

- 30,000 Lux @ 0.5m
- continuously dimmable illumination field with head connection vertical/horizontal
- colour temperature fixed 4400 K



#### Optimal vision

thanks to homogeneous, spot-like illumination field

**medicinelaser**  
LATAM

[www.medicinelaser.com](http://www.medicinelaser.com)



ARGENTINA  
Av. Doctor Honorio  
Pueyrredón 707, Caballito,  
CABA, Argentina.  
M: +54 11 3404 7158

CHILE  
Sucre 1355, Ñuñoa,  
Santiago, Chile.  
O: +56 2 2269 1981  
M: +56 9 4409 7105

ECUADOR  
Av. República 1058,  
Edificio La Fontana,  
Quito, Ecuador.  
M: +593 99 120 4597

PERÚ  
Calle Almirante Lord  
Cochrane 322,  
Miraflores, Lima, Perú.  
O: + 51 1 704 80 14  
M: +51 983 552 813

+ OFICINAS CENTRALES

[ventas@medicinelaser.com](mailto:ventas@medicinelaser.com)