

# OPTICLUX Hand



medicinelaser<sup>+</sup>  
LATAM



## Strikingly different. OPTICLUX Hand

The OPTICLUX Hand magnifier luminaire with ergonomic handle provides support for demanding visual tasks. This is ensured by the large dimmable light field in daylight quality, excellent colour rendering properties and long battery life. Dermatologists and veterinarians can look forward to the powerful Wood's Light, which is available as an option in combination with the daylight mode.



OPTICLUX Hand	technical data
<b>luminous power *</b>	8,000 Lux @ 0.15 m
<b>colour temperature *</b>	<ul style="list-style-type: none"> <li>• 6500 K</li> <li>• optional: 365 nm (Wood's Light)</li> </ul>
<b>lens dimension</b>	Ø 16 cm
<b>colour rendering index</b>	Ra > 95
<b>usage</b>	touch panel on handle
<b>class of protection</b>	II
<b>protection rating</b>	IP 20



special features
<ul style="list-style-type: none"> <li>• dimmable</li> <li>• rechargeable battery with USB-C charging interface</li> <li>• Wood's Light mode (optional)</li> <li>• lens cover</li> <li>• 1.9x magnification</li> </ul>
<b>OPTICLUX Hand</b>
<b>variations</b>
battery operated handheld version

# OPTICLUX Hand – variations



## OPTICLUX Hand

Daylight & Wood's Light



## OPTICLUX Hand + Table charger

Setprice available



### First class light

The OPTICLUX hand magnifier light delivers unique light output of 8,000 lux at a distance of 15 cm. Excellent colour rendering properties > Ra 95 ensure ideal vision.

### Large field of view

The scratch-resistant plastic lens with 160 mm diameter and 3.5 diopters is ideal for full-body examinations and for assessing dermatological diseases on patients.

## Product highlights



### Optional with Wood's Light

two lighting modes in one device: optional Wood's Light for examination of fungal or bacterial infection on the skin or scalp and pigmentary changes

For detailed description of Wood's Light see page 110.



### Ergonomic

optimal handling for right- and left-handed



**Table charging station**  
for storage and charging, white powder-coated with a sturdy aluminium base



**Easy to operate**  
touch panel on the handle with comfortable, stepless dimming



**Battery operated**  
up to 3 hours battery runtime, flexible charging through USB-C interface

**medicinelaser**  
LATAM  
[www.medicinelaser.com](http://www.medicinelaser.com)



ARGENTINA  
Av. Doctor Honorio  
Pueyrredón 707, Caballito,  
CABA, Argentina.  
M: +54 11 3404 7158

CHILE  
Sucre 1355, Ñuñoa,  
Santiago, Chile.  
O: +56 2 2269 1981  
M: +56 9 4409 7105

ECUADOR  
Av. República 1058,  
Edificio La Fontana,  
Quito, Ecuador.  
M: +593 99 120 4597

PERÚ  
Calle Almirante Lord  
Cochrane 322,  
Miraflores, Lima, Perú.  
O: + 51 1 704 80 14  
M: +51 983 552 813

[ventas@medicinelaser.com](mailto:ventas@medicinelaser.com)